

Brain Computer Interfaces: Fundamentals and Application

Lecture Week 13. ERP-based BCIs

Chun-Hsiang (Michael) Chuang



Today's Agenda

ERP-based BCIs

- Stimulus-evoked activity
- P300
- SSVEP
- RSVP

Online course

<https://meet.google.com/fba-uxbj-md>

EEG-based BCIs

Modulating the response of large neural populations

- Subject training over a period of time → *self-paced (asynchronous BCI)*
voluntarily initiate control
utilize some form of imagery (motor or cognitive)
- External stimuli → *Stimulus-based (synchronous BCIs)*
stereotypical brain responses

Which is easier to be built?

Two major types of BCIs

Invasive BCIs


Noninvasive BCIs

- Only record from the brain and translate the neural data into control signals for output devices.
- Only stimulate the brain and cause certain desired patterns of neural activity in the brain.
- Both record and stimulate the brain

Oscillatory Potentials and ERD

Performs movement or imagines performing a movement

- mu (8–12 Hz) or beta band (13–30 Hz) decreases
i.e., event-related desynchronization (ERD)



Brain Responses Useful for Building BCIs

Conditioned Responses

Population Activity

Imagined Motor and Cognitive Activity

Stimulus-Evoked Activity

- Brain activity generated in response to certain stimuli
- P300 (or P3) in EEG, an event-related potential (ERP), over the parietal lobe
- The exact neural mechanisms responsible for the P300 are as yet unclear

Other common types of evoked potentials

Steady state visually evoked potential (SSVEP)

- Visual cortex
- A visual stimulus flickering at a certain frequency
- Peaks in the power spectrum at the same frequency

N100 (or N1)

- Frontal and Central regions
- Usually followed by P200 (N1-P2 complex)
- Sudden and loud noise

Other common types of evoked potentials

N400

- Central and Parietal sites
- Incongruent inputs
- Similar to an *error potential* (ErrP)

Stimulus-Evoked Potentials

A “rare” stimulus can evoke a potentials with a positive peak at 300 ms after the stimulus is presented.

- auditory, visual, or somatosensory forms.

visually evoked potentials (VEPs)

- SSVEP

auditory evoked potentials (AEPs)



somatosensory evoked potentials (SSEPs)

The P300 Potential (250 to 500 ms)

Stimulus must be rare and unpredictable but relevant.

Amplitude of P300

- how relevant the stimulus is
- varies inversely with the probability of stimulus
- decision making




P3a (or novelty P3)

maximum amplitude over
frontal/central electrode sites

peak latency in the range of
250-280 ms

associated with
engagement
of attention (especially
the **orienting**, involuntary
shifts to changes in the
environment)

processing of **novelty**




P3b

peaks at around 250–500 ms,
over **parietal** brain areas.

The P3b has been a prominent
tool used to study cognitive
processes, especially psychology
research on **information
processing**.

The P3b can also be used to
measure how demanding a task
is on cognitive workload.




BCI “speller based on “Oddball paradigm”

Farwell and Donchin (1988)

- 26 letters + symbols/commands
- Focusing attention on the target letter
- Rows and columns of the matrix repeatedly flash in random order

Choose one letter or command

| | | | | | |
|---|---|---|---|-----|------|
| A | G | M | S | Y | * |
| B | H | N | T | Z | * |
| C | I | O | U | * | TALK |
| D | J | P | V | FLN | SPAC |
| E | K | Q | W | * | BKSP |
| F | L | R | X | SPL | QUIT |




BCI “speller based on “Oddball paradigm”

Containing a subject's chosen letter → a large P300 generated.

Linear discriminate analysis classifier

Speed-accuracy trade-off:

Higher the number of flashes, the better the accuracy, but ...



Steady State Visually Evoked Potential (SSVEP)

continuously fluctuating stimuli (with repetition rate > 5 Hz)



Focus attention on the button corresponding to your choice.

EEG in the occipital region oscillates at the stimulus frequency.


So, how to extract EEG features?

Steady State Visually Evoked Potential (SSVEP)

SSVEP evoked by 7 Hz visual stimulation

FFT

Peaks at 7 Hz and its harmonics





Steady State Visually Evoked Potential (SSVEP)

6-14 Hz

FFT was performed every 0.3s to compute the amplitude spectrum


Sum(amplitude and harmonics) vs mean of amplitude spectrum

SSVEP-based method can obtain the highest information transfer rate



BCI using SSVEP

[4-command SSVEP BCI demo](#) (by paweu)



BCI using SSVEP

SSVEP-Based Brain-Computer Interface (BCI) System for Real-Time Application (by University of Malta)




Exercise: SSVEP

Step 1. Select one of SSVEP datasets

<https://www.mamem.eu/results/datasets/>

Step 2. Apply FFT to analyze the power spectrum of EEG activity in response to the visual stimuli

Step 3. Try to build a classifier that can help you distinguish five different frequencies



MAMEM
Multimedia Authoring & Management
using your Eyes & Mind

Toggle navigation through eyes & mind

Home The Project Results Use Cases News

Datasets of the project

Available datasets related to the MAMEM Project.


| Title | Description | Download |
|----------------------|--|--------------------------|
| EEG SSVEP Dataset I | EEG signals with 256 channels captured from 11 subjects executing a SSVEP-based experimental protocol. Five different frequencies (6.66, 7.50, 8.57, 10.00 and 12.00 Hz) presented in isolation have been used for the visual stimulation. The EGI 300 Geodesic EEG System (GES 300), using a 256-channel HydroCel Geodesic Sensor Net (HCGSN) and a sampling rate of 250 Hz has been used for capturing the signals. | Download |
| EEG SSVEP Dataset II | EEG signals with 256 channels captured from 11 subjects executing a SSVEP-based experimental protocol. Five different frequencies (6.66, 7.50, 8.57, 10.00 and 12.00 Hz) presented simultaneously have been used for the visual stimulation. The EGI 300 Geodesic EEG System (GES 300), using a 256-channel HydroCel Geodesic Sensor Net (HCGSN) and a sampling rate of 250 Hz has been used for capturing the signals. | Download |

Rapid Serial Visual Presentation (RSVP)

Examine the temporal characteristics of **attention**

The RSVP paradigm requires participants to look at **a continuous presentation of visual items** which is around 10 items per second.

The **targets** are placed inside this stream of continuous items. They are separate from the rest of the items known as **distractors**.



Example

[RSVP Keyboard TM](#) (by CogSysLab atNEU)




A Benchmark Dataset for RSVP-Based Brain–Computer Interfaces

Zhang et al. (2020), <http://bci.med.tsinghua.edu.cn/download.html>

Images: Computer Science and Artificial Intelligence Library of MIT University

target (with human) vs. non-target (without human)






Auditory Evoked Potentials

Replace visual stimuli with auditory stimuli → auditory oddball paradigm



- 50 ms square-wave beeps
- Frequent non-target beeps + occasional target beeps
- Paying attention to the target stimuli (by counting)
- 39-ch EEG signals, unmixed using ICA, and classified using a linear SVM → 5-15% error rate

Tutorial: ERP





■ Visualize data measures



■ Manually plot a topography


```
chanlocs = struct('labels', {'FP1' 'FP2' 'F7' 'F3' 'FZ' 'F4' 'F8' 'FT7' 'FC3'  
'FCZ' 'FC4' 'FT8' 'T3' 'C3' 'CZ' 'C4' 'T4' 'TP7' 'CP3' 'CPZ' 'CP4' 'TP8' 'T5'  
'P3' 'PZ' 'P4' 'T6' 'O1' 'OZ' 'O2'});
```

```
chanlocs = pop_chanedit(chanlocs);
```



```
topoplot(EEG.data(:, 1000), chanlocs);
```

■ Visualize data measures

headplot.m





Plot data spectrum and maps




Plot channel ERPimage

ERP: event-related potentials




Reaction Time

EEG.event.latency






Plot a RT-sorted ERPimage









Visualize data measures

36



Visualize data measures





Co-register channel locations with head mesh and compute a mesh spline file (each scalp montage needs a headplot() spline file)

- Use the following spline file or structure
 Or (re)compute a new spline file named:
 3-D head mesh file
 Mesh associated channel file
 Talairach-model transformation matrix

| | | | | |
|-----------|----------|---------|-----------|---------|
| -0.355789 | -6.33688 | 12.3705 | 0.0533239 | 0.01874 |
|-----------|----------|---------|-----------|---------|

Plot interpolated activity onto 3-D head

Making headplots for these latencies (from -1000 to 6996 ms):


Plot title:

Plot geometry (rows,columns): (Default [] = near square)



Other headplot options (See >> help headplot):

ERP scalp maps of dataset:s05_061019m res


Visualize data measures



2D vs 3D



Possible Solution – Independent Component Analysis



The signal recorded can be considered as a linear mixture of underlying neural generators

$$A \times S(t) = X(t)$$

Assume that Sources are statistically independent from each.


$$W \times X(t) = \hat{S}(t)$$

**Unmixing
Matrix**

Independent Component Analysis finds an unmixing matrix W to maximize some objective function that measures **independence of the estimated components**

Independent Component Analysis (ICA)

Unmixes the uncertainty of **WHAT** from that of **WHERE!**



$$\begin{aligned} \mathbf{W} \times \mathbf{X}(t) &= \hat{\mathbf{s}}(t) \\ \mathbf{X}(t) &= \mathbf{W}^{-1} \times \hat{\mathbf{s}}(t) \\ &\downarrow \\ &\text{Inverse of} \\ &\text{Unmixing Matrix} \end{aligned}$$

Given n sensors,

$$\begin{aligned} \mathbf{W}^{-1} &= \Lambda \\ &= \begin{bmatrix} \lambda_{11} & \cdots & \lambda_{1i} & \cdots & \lambda_{1n} \\ \lambda_{21} & \cdots & \lambda_{2i} & \cdots & \lambda_{2n} \\ \vdots & \ddots & \vdots & \ddots & \vdots \\ \lambda_{n1} & \cdots & \lambda_{ni} & \cdots & \lambda_{nn} \end{bmatrix} \end{aligned}$$


Independent Component Analysis (ICA)

$$\mathbf{X}(t) = \mathbf{W}^{-1} \times \hat{\mathbf{S}}(t)$$


$$\mathbf{W}^{-1} = \Lambda$$

$$= \begin{bmatrix} \lambda_{11} & \cdots & \lambda_{1i} & \cdots & \lambda_{1n} \\ \lambda_{21} & \cdots & \lambda_{2i} & \cdots & \lambda_{2n} \\ \vdots & \vdots & \vdots & \vdots & \vdots \\ \lambda_{n1} & \cdots & \lambda_{ni} & \cdots & \lambda_{nn} \end{bmatrix}$$

$$\begin{bmatrix} X_1(t) \\ \vdots \\ X_n(t) \end{bmatrix} = \mathbf{W}^{-1} \times \begin{bmatrix} \hat{\mathbf{S}}_1(t) \\ \vdots \\ \hat{\mathbf{S}}_n(t) \end{bmatrix}$$



EEGLAB



Kurtosis will be calculated initially every 1 blocks using 6000 data points.

Decomposing 1517 frames per ICA weight ($(900)^2 = 1366000$ weights), Initial learning rate will be 0.001, block Learning rate will be multiplied by 0.98 whenever angledelta ≥ 60 deg.

More than 32 channels: default stopping weight change 1E-7

Training will end when wchange $< 1e-06$ or after 512 steps.

Online bias adjustment will be used.

Removing mean



Final training


Computing the

Starting weight

Sphering the

Beginning ICA

Warning: Usin



generators. This syntax is not recommended. See [Replace Discouraged Syntaxes of rand and randn to use RNG](#)

to replace the old syntax.

> In [runica](#) (line 830)

 In [pop_runica](#) (line 399)

Lowering learning rate to 0.0009 and starting again.

Lowering learning rate to 0.00081 and starting again.


Lowering learning rate to 0.000729 and starting again.

```
ommand Window
step 48 - lrate 0.000019, wchange 0.60102712, angledelta 137.4 deg
step 49 - lrate 0.000019, wchange 0.45452005, angledelta 169.7 deg
step 50 - lrate 0.000018, wchange 0.48941163, angledelta 170.1 deg
step 51 - lrate 0.000018, wchange 0.04110124, angledelta 75.0 deg
step 52 - lrate 0.000018, wchange 0.40831810, angledelta 115.8 deg
step 53 - lrate 0.000017, wchange 0.73758322, angledelta 167.0 deg
step 54 - lrate 0.000017, wchange 0.05327091, angledelta 144.9 deg
step 55 - lrate 0.000016, wchange 0.02706516, angledelta 124.5 deg
step 56 - lrate 0.000016, wchange 0.01631103, angledelta 104.6 deg
step 57 - lrate 0.000016, wchange 0.02201465, angledelta 120.1 deg
step 58 - lrate 0.000016, wchange 0.33487233, angledelta 113.6 deg
step 59 - lrate 0.000015, wchange 0.10354193, angledelta 33.3 deg
step 60 - lrate 0.000015, wchange 0.01142147, angledelta 95.9 deg
step 61 - lrate 0.000015, wchange 0.01723309, angledelta 100.1 deg
step 62 - lrate 0.000015, wchange 0.01847892, angledelta 145.3 deg
step 63 - lrate 0.000014, wchange 0.56760624, angledelta 94.6 deg
step 64 - lrate 0.000014, wchange 0.03602543, angledelta 110.7 deg
step 65 - lrate 0.000014, wchange 0.03381547, angledelta 142.6 deg
step 66 - lrate 0.000013, wchange 0.01898318, angledelta 114.5 deg
step 67 - lrate 0.000013, wchange 0.00898991, angledelta 122.3 deg
step 68 - lrate 0.000013, wchange 0.01130301, angledelta 123.4 deg
x>> |
```

EEG.icaact?

| | |
|-------------|--------------------------|
| xmin | -1 |
| xmax | 1.9960 |
| times | <i>1x750 double</i> |
| data | <i>30x750x683 single</i> |
| icaact | [] |
| icawinv | <i>30x30 double</i> |
| icasphere | <i>30x30 double</i> |
| icaweights | <i>30x30 double</i> |
| icachansind | <i>1x30 double</i> |





Unmixing Matrix

$$\mathbf{A} \times \mathbf{S}(t) = \mathbf{X}(t)$$

$$\mathbf{W} \times \mathbf{X}(t) = \hat{\mathbf{S}}(t)$$

$$\mathbf{X}(t) = \mathbf{W}^{-1} \times \hat{\mathbf{S}}(t)$$

Inverse of
Unmixing Matrix

Given n sensors,

$$\mathbf{W}^{-1} = \Lambda$$

$$= \begin{bmatrix} \lambda_{11} & \cdots & \lambda_{1i} & \cdots & \lambda_{1n} \\ \lambda_{21} & \cdots & \lambda_{2i} & \cdots & \lambda_{2n} \\ \vdots & \ddots & \vdots & \ddots & \vdots \\ \lambda_{n1} & \cdots & \lambda_{ni} & \cdots & \lambda_{nn} \end{bmatrix}$$

| | |
|-------------|--------------------------|
| xmax | -1 |
| times | 1.9960 |
| data | <i>1x750 double</i> |
| icaact | <i>30x750x683 single</i> |
| icawinv | [] |
| icasphere | <i>30x30 double</i> |
| icaweights | <i>30x30 double</i> |
| icachansind | <i>1x30 double</i> |

icaact??

The Infomax ICA algorithm first "sphere" the data (remove correlation) and then performs ICA.

So the inverse weight matrix EEG.icawinv is equal to

$$\text{EEG.icawinv} = \text{pinv}(\text{EEG.icaweights} * \text{EEG.icasphere});$$

$$\mathbf{W} \times \mathbf{X}(t) = \hat{\mathbf{S}}(t)$$

$$\mathbf{X}(t) = \mathbf{W}^{-1} \times \hat{\mathbf{S}}(t)$$

$$\mathbf{X}(t) = \mathbf{W}^{-1} \hat{\mathbf{S}}(t) \quad \mathbf{W} \times \mathbf{X}(t) = \hat{\mathbf{S}}(t)$$

**Inverse of
Unmixing Matrix**



Channels domain

vs.



Component domain



Frontal



Left motor



Parietal

